



# LambGro (22/28)

Milk Replacer for Rearing Replacement Ewe Lambs & Meat Lambs  
[www.GroberNutrition.com](http://www.GroberNutrition.com)

## Product Overview

Grober LambGro milk replacer is designed for replacement feeding of multiple birth lambs, as well as for feeding orphaned lambs. The optimal concentration matches the nutritional profile of ewe's milk.

Features	Benefits
100% milk protein source	Protein and amino acid balance for optimal growth and metabolic controls
Skim milk powder	Provides highly digestible nutrients for growth and health
Optimal Vitamin E & C	Antioxidants to promote health and immune function
Machine feeding	Can be used in automatic lamb feeding machines

## Lamb Management Tips

- ✓ Place lambs in a clean and dry area immediately after birth
- ✓ Feed LambGro/KidGro™ Colostrum as soon as possible after birth in 3 feedings within 18 hours of life
  - ✓ Feed by bottle or use a feeding tube if necessary
- ✓ Ensure the lamb pen is free of draft and is clean and warm
- ✓ Start feeding milk replacer daily up until weaning
- ✓ Free-choice water should be available at all times
- ✓ Multiple feedings/day recommended to reduce incidence of bloat

## Shelf Life and Storage

Store in a cool, dry place up to 1 year from date of manufacture. Store opened bags in a sealable container.

## Mixing Directions

1. Weigh 190g of powder per litre of solution using a scale.
2. Mix with 0.5 litre hot water (60°C) for 3 minutes.
3. Add 0.45 litre cool water and mix for 1 minute (total 1 litre).
4. Feed the milk solution at a temperature of 38-40°C (body temperature) (thermometer provides best results).
5. Clean all feeding utensils and equipment after each use.

## Feeding Schedule

Age of Lamb	Amount Per Feeding (ml)	Number of Feedings/Day
Day 1 (colostrum)	200	4
The newborn lamb should receive 150mL/kg bodyweight of colostrum (approx. 600mL), within six hours of birth. (*LambGro/KidGro™ Colostrum Replacer)		
Day 2-4	200	4
Day 5-7	300	4
Day 8-14	500	3
Week 3-4	600	3
Week 5-6	700	2
Week 7-8	700	1



Package Size

20KG 10KG 5KG